

AKED.bzh



FINGERFIX PLIERS


REF.: 0227 02 00

FINGERFIX Pliers

REF.: 0227 02 00



Length: 16.3 cm - Weight: 197 g



FINGERFIX Pliers can bend and cut K-wires measuring 1 mm to 1.6 mm in diameter. FINGERFIX Pliers were designed in order to bend a K-wire without applying mechanical stress while traversing bone.

BENDING

Positioning:

The K-wire is inserted laterally in one of the slots of the lower jaw.


Bending :



Bending a 1.5 mm \varnothing Kirschner wire.

During the bending operation, the pin is pushed backwards by the nose of the upper jaw and guided by a small groove.

The right lateral slot of the plier's lower jaw can hold pins measuring 1.6 mm or 1.2 mm in diameter. The left hand slot, which is smaller, can hold pins measuring 1 mm or 1.2 mm in diameter.



\varnothing : 1 et 1,2mm

crimp

\varnothing : 1,2 et 1,6mm

N.B.

- » The pin need only protrude 20 mm from the surface of the skin to be bent effortlessly by the FINGERFIX Pliers.
- » To fully apprehend the fact that bending the pin with FINGERFIX Pliers does not affect it, hold the pin between thumb and index finger placed in contact with the pliers.

CUTTING



The pin is fully inserted into the cutting aperture (like a cigar cutter).

FINGERFIX CUT



STANDARD CUT



FINGERFIX Pliers cuts by shearing (like a cigar cutter) and not by crushing.

As a result, the cutting surface on the end of the pin is clean and burr-free, unlike the surface with standard cutting instruments.

Photographs of cuts on a 1.6 mm diameter Kirschner wire.

- » Guaranteed for 5 years
- » 1 year for the jaws
- » Stainless steel.
- » Copyright ©