# KIROFIX

*Main malléable / malleable hand /  
Formbare Hand / apoyamano maleable /  
mano malleabile / mão maleável /  
vormbare hand*

## KIROFIX BALLAST

*lest pour main malléable / ballast for malleable  
hand / Ballast für formbare Hand / lastre para el  
apoyamano maleable / zavorra per mano  
malleabile / estabilizador para mão maleável /  
ballast voor vormbare hand*

## KIROFIX RING

*Élastique de fixation de doigt pour KIROFIX  
Elastic finger fixation for KIROFIX  
Gummi zur Fingerfixierung für KIROFIX  
Elástico de fijación de dedos para KIROFIX  
Elastico di fissaggio del dito per KIROFIX  
Elástico de fixação de dedo para a KIROFIX  
Elastische vingerfixatie voor KIROFIX*

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# EN

## DEVICE AND INDICATION FOR USE

- **KIROFIX [KIROFIX Large, KIROFIX Small and KIROFIX Medium]** | KIROFIX, a reusable malleable silicone hand, is designed to hold the patient's hand in the desired position during an orthopedic surgery. KIROFIX Large is intended to be used on adult size hands, KIROFIX Small is intended to be used on children size hands, KIROFIX Medium can be used on all size of hands. KIROFIX is compatible with MRI (Magnetic Resonance Imaging).

- **KIROFIX BALLAST** | KIROFIX Ballast is a reusable stainless-steel accessory which allows to ballast the assembly in order to stabilise KIROFIX in the desired position.

- **KIROFIX RING** | KIROFIX RING is a reusable silicon component of KIROFIX Medium designed to be used as an elastic finger fixation. KIROFIX RING can be also used as an accessory with all KIROFIX range for the same function.

These devices have been designed to facilitate hand surgery. These devices which are supplied non-sterile must undergo an automatic cleaning step followed by a steam sterilization step before use.

## WARNINGS

- Do not use the instrument when it comes out of the sterilization cycle and is still warm : there is a risk of premature wear.
- KIROFIX must be completely unfolded and well flattened before cleaning and sterilization otherwise its lifetime will be divided by 2.
- Do not let humidity remain on the KIROFIX Ballast's laser marking. Otherwise, if a brownish spot appears, rub it away vigorously with a soft cloth.
- Do not store near corrosive products (ex: chlorine).
- Do not use the instruments if they are damaged. Example: If a defect is found in the KIROFIX's plugs, it is non-compliant and no longer usable.
- Do not clamp KIROFIX RING with a metal instrument, this can make small cuts in the silicone that may cause failure later on.
- Failure to follow good practices for the cleaning and sterilization of surgical instruments may lead to infections or device malfunction.
- Do not use KIROFIX and its accessories for any use other than those defined by the manufacturer.

## LIMITATIONS ON REPROCESSING

The instruments can be reused, however, their lifespan is limited. The lifespan of KIROFIX and its accessories depends on the individual use (wear, damage). The devices can be reused if the conditions of cleaning and sterilization provided by ASKORN are followed. The mechanical strength of these devices has been confirmed by ASKORN after the equivalent of 225 cleaning/sterilization cycles.

In order to avoid damaging of KIROFIX and its accessories, it is forbidden to use solutions containing: chlorine, soda and aromatic derivatives compounds (Ether, Chloroform, Benzene, Xylene etc.) or any other abrasive product when the cleaning and the decontamination these devices. In the case of suspected contamination by non-conventional

transmissible agents (NTCA: prion disease), the devices cannot be reused and must be destroyed according to current regulations.

## **INSTRUCTIONS**

### **1. Putting the new product into circulation (out of the packaging):**

Remove the devices from their packaging.

It is recommended to clean the new device before the first sterilisation according to point 5.

### **2. Initial processing at point of use (after surgical use):**

Cleaning must be carried out immediately after each use of the medical device to prevent corrosion and dirt adhesion. Use a cloth to remove any excess dirt.

### **3. Containment and transport**

The mode of transport must be secured, no additional requirements.

### **4. Preparation before automatic cleaning**

It is recommended to clean new instruments before their initial sterilization and immediately after each use in order to avoid dirt adhesion.

- Unfold KIROFIX
- Separate the KIROFIX Ballast from KIROFIX
- Separate the strap from the KIROFIX Medium
- Separate KIROFIX RING

KIROFIX, KIROFIX Ballast and KIROFIX RING will be cleaned like surgical instruments.

### **5. Automated cleaning**

The used detergent is an enzyme-based or alkaline based.

Use a cleaning program similar to this one described below:

1. Pre-cleaning  
Tap water / Cold water / 2 min.
2. Drying
3. Cleaning  
Tap water + Medizyme enzymatic ou Mediclean forte, alkaline (Neodisher\*) 0.5% / 55°C / 5 min.
4. Drying
5. Intermediate rinsing  
Demineralized water / Cold water / 3 min.
6. Drying
7. Final rinse  
Demineralized water / Cold water / 2 min.
8. Drying

\* Follow the detergent manufacturer's instruction. Note: A high temperature disinfection step can take place after the final rinse to prevent the fixation of any remaining proteins.

It is recommended to use a washer/disinfector conforms to ISO 15883.

2 automatic cleaning cycles can be performed prior to sterilization if necessary. After cleaning, always check that the devices are dry, clean, in good working order and absence of deterioration.

### **6. Automated drying**

Dry KIROFIX, KIROFIX RING and KIROFIX Ballast according to the drying cycle of the washer. If necessary, additional manual drying can be performed using lint free towel.

## **7. Maintenance**

Non-applicable.

## **8. Controls and functional tests**

After each cleaning and drying, check the dryness, cleanliness, proper functioning and absence of deterioration of the devices.

Before each use, carry out the following tests in order to verify the conformity of the device:

- **KIROFIX**

Visual inspection of the device to identify any damage, wear and/or signs of degradation: hole, crack, air pocket in the hand, defect on KIROFIX's plugs, broken finger under the silicon cover.

- **KIROFIX BALLAST**

Visual inspection of the device to identify any damage wear and/or signs of degradation : corrosion

- **KIROFIX RING**

Visual inspection of the device to identify any damage wear and/or signs of degradation: crack

## **9. Packaging**

Standard individual packaging compatible with steam sterilization method must be used (i.e. sealed bags). Make sure the pouch is enough large to hold the unfolded instrument without placing stress on the pouch's seal. The KIROFIX Ballast is always packaged separately.

## **10. Sterilization**

KIROFIX must be unfold and flattened before sterilization. KIROFIX Ballast, KIROFIX et RING are sterilized in a vacuum autoclave according to one of the following protocols:

### **Protocol A**

- 3 pre-vacuum phases
- Sterilization temperature: 132°C
- Holding time: 4 minutes
- Drying time: 20 minutes

### **Protocol B**

- 3 pre-vacuum phases
- Sterilization temperature: 134°C
- Holding time: 3 minutes
- Drying time: 20 minutes

The instructions of the autoclave's manufacturer must be followed.

Ensure that the surfaces of the device are not covered during sterilisation.

Only the cleaning and sterilization procedures defined in this IFU have been validated. However, the healthcare facility can use its standard disinfection and sterilization protocol, as long as this protocol has been validated for the disinfection and sterilization of surgical instruments (e.g. sterilization in instrument holders, manual cleaning, other cleaning and sterilization parameters).

## **STORAGE CONDITIONS**

KIROFIX, KIROFIX Ballast and KIROFIX RING must be stored in a dry and clean and dry environment, at room temperature. Do not store the instruments near products that may have a corrosive action (ex: chlorine for KIROFIX Ballast).

## **DISPOSAL**

Any instrument that is broken, damaged or at the end of its life must not be used. It must be destroyed like any medical waste according to current regulations.

## **ADDITIONAL INSTRUCTIONS**

### **1. General Information**

KIROFIX, KIROFIX Ballast and RING are provided non-sterile. Devices must be used by surgeon and qualified personnel. Instruments must be used only for the surgical purpose and type for which they have been designed (see indications). Any failure to respect these rules cancels the guarantee:

- KIROFIX (KIROFIX Large, KIROFIX Small, KIROFIX Medium): 1 year
- KIROFIX Ballast: 5 years
- RING: 1 year

If for any reason, some instrument needs to be returned to us, it must be first decontaminated. ASKORN reserves the right to request a decontamination proof.

### **2. Precautions for use**

- KIROFIX, KIROFIX Ballast and RING must be used in accordance with this Instruction for Use.
- Read all the sections of this IFU before use.
- The instruments should be cleaned immediately after use.
- Before use, inspect devices to verify its condition and proper functioning.

### **3. Incident declaration**

Any serious incident that has occurred in relation to KIROFIX and its accessories should be reported to ASKORN, and the competent authority of the Member State in which the user and/or patient is established.

Reminder: A serious incident or risk of serious incident is any incident or risk of incident involving a medical device which led or might have led to death or serious deterioration in state of health of patient, end-user and third party.

*These instructions have been validated by the medical device's manufacturer as being appropriate for preparing a medical device for reuse. The department is still responsible for ensuring the sterilization process achieves the intended result by involving the sterilization facility's staff, materials and equipment. This generally requires routine verification and control of the procedures. Similarly, any deviation from the instructions provided should be evaluated to determine the effectiveness and possible adverse consequences.*

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