FINGERFIX NON STERILE

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INSTRUCTION FOR USE

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MANUAL DE INSTRUCCIONES

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IFU221124 Rev3

DESCRIPTION

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DESCRIPTION

The FINGERFIX NON STERILE is a digital dynamic external fixator. It consists of 2 plastic connector rods and 2 stainless steel springs. In order to guide the double spring system, it is necessary to use 2 pins of Ø1.2 mm and the pin of Ø1.5 mm in diameter. This third U-shaped pin connects the two distraction mechanisms.

To ensure optimal protection, the FINGERFIX NON STERILE is also composed of 4 pin protective tips and a single-use reward guide to ensure the parallelism of the two 1.2 mm pins. For ideal use, the NON STERILE FINGERFIX is combined with the FINGERFIX PLIERS which allow the pins to be pinched at 90° (without mechanical stress on the bone) and to cut them cleanly as well as to ensure the crimping of the protective pin tips.

The device is designed for a maximum duration of 60 days.

INDICATIONS

The FINGERFIX NON STERILE is a dynamic external fixator designed to treat dislocations of complex, stable and unstable dislocations, fracture luxations and pilon fractures at the level of the proximal interphalangeal joint of the fingers of the hand. It is designed to treat contractures of the interphalangeal joints with associated loss of mobility. The device restores finger alignment and mobility while allowing early mobilization and anatomically normal movements.

For children, the surgeon will determine whether FINGERFIX NON STERILE can be used, taking into account age, morphology and ability to follow post-operative follow-up instructions.

CONTRAINDICATIONS

- Chronic or acute local or systemic infection, fever.
- Known allergy to chromium or any other component of stainless steel.
- Insufficient tissue coverage.
- Infectious disease, metabolic or systemic disorder.
- Drug or alcohol addiction.
- Severe osteoporosis, local bone tumour.
- Poor nutritional supply of tissues, insufficient bone stock.
- Obvious intellectual incapacity of the patient to follow the surgeon's recommendations: non-cooperative patient or patient with a mental illness.
- Intense physical activity.

CAUTION

FINGERFIX NON STERILE cannot be used with MRI because it may cause skin burns.

It is advisable to clearly explain to the patient the objectives and the operating mode of the distraction so that the treatment proceeds for the best. Tell him the warning signs that will make him contact his surgeon so that corrective action can be implemented: loosening, disappearance of traction, breakage, unusual pain, infection, purulent discharge.

CLEANING

Automated cleaning

Steps :

1. Pre-cleaning

Tap water: 16°C ± 2°C — 3 min

2. Drain

3. Cleaning: Tap water: 60°C — 5 min

Products recommended by ASKORN MEDICAL: Deconex prozyme alka (Borer). Concentration: 0.5%

4. Drain

5. Intermediate rinse: Deionized water: 45°C — 1 min

6. Drain

7. Final rinse: Deionized water: 20°C ± 2°C - 0.5 ml/l — 1 min

8. Drain

9. Thermal disinfection: Demineralised water: 90°C — 5 min

10. Drying: No water used: 90°C — 15 min

2 automatic cleaning cycles can be performed before sterilization when necessary.

An alkaline detergent (pH up to 11) can also be used (e.g. Neodisher Mediclean Forte). In this case, follow the detergent manufacturer's recommendations.

After cleaning, always check the dry condition, cleanliness, correct operation and absence of deterioration of the device(s).

We recommend using a washer-disinfector compliant with the ISO 15883 series of standards.

Manual cleaning

Steps :

1. Pre-cleaning

Tap water: $16^{\circ}\text{C} \pm 2^{\circ}\text{C}$ — 5 min

2. Pre-cleaning

Ultrasonic bath: 20°C — 5 min

Products recommended by ASKORN MEDICAL: Deconex

Prozyme Active (Borer). Concentration: 0.2%

3. Rinse

Tap water: $16^{\circ}\text{C} \pm 2^{\circ}\text{C}$ — 20 sec

Use a water jet gun with approximately 2 bar static pressure.

2 manual cleaning cycles can be performed before sterilization when necessary.

After cleaning, always check the dry condition, cleanliness, correct operation and absence of deterioration of the device(s).

Drying

• Automatic drying

Dry the items according to the washer/disinfector drying cycle as detailed above. If needed, additional manual drying can be done with a lint-free towel.

• Manual drying

Dry with a lint-free towel.

PACKAGING

Standard individual packaging compatible with steam sterilization must be used (i.e. sealed pouches, box with filter and valve or crepe paper). It must be ensured that the bag is large enough to contain the devices without stressing the seal of the bag.

STERILIZATION

Place in sterilizer. Items should be sterilized by vacuum autoclave according to one of the following protocols:

Protocol A

- 3 pre-vacuum phases.
- Sterilization temperature of 132°C .
- Duration of temperature maintenance: 4 min
- Drying time: 18 min

Protocol B

- 3 pre-vacuum phases.
- Sterilization temperature of 134°C .
- Duration of temperature maintenance¹: 3 min
- Drying time: 18 min

¹Items are compatible with an 18-minute exposure to 134°C .

The autoclave manufacturer's instructions should be followed. Ensure that the inner and outer surfaces of the device are freely accessible during sterilization. Only the cleaning and sterilization processes defined in these instructions for use have been validated by the manufacturer, however the care unit concerned may use its standard disinfection and sterilization protocol, provided that said protocol is validated under its responsibility.

If the device is not used before the sterile state expires, it is possible to reprocess the device according to the instructions in this manual up to a total of 5 cleaning and sterilization cycles.

SURGICAL PROCEDURE

The implementation of appropriate surgical procedures and techniques is the sole responsibility of the medical profession.

PATIENT INFORMATION

Any surgical procedure has risks and can be accompanied by complications. It is up to the surgeon to clearly explain the benefit/risk ratio of the proposed technique. Possible complications and warnings should be discussed with the patient before the decision to proceed with surgery is made. The possibility of interrupting treatment at any time should also be discussed with the patient.

POSSIBLE COMPLICATIONS AND WARNINGS

Complications that may be associated with all surgical procedures should be discussed with the patient (infection, hematoma, fluid accumulation, loss of sensation, pain, stiffness, deformity-deviation of the finger, algodystrophy, subluxations, possible allergic reaction to drugs used during intervention, nerve damage, intolerance of the patient to a foreign body, and poor healing).

DISPOSABLE

The devices are intended for single use only.

Explanted products should not be reused because re-cleaning and re-sterilization procedures may not safely remove all biological residues such as blood, tissues and other materials that may contain resistant pathogens.

PACKAGING

The integrity of the device is only guaranteed if the packaging is intact. Peel-off labels to ensure product traceability are provided in the packaging.

HOW TO USE

Implementing correct surgical procedures and techniques is a matter for the medical profession. Each surgeon should assess the suitability of the procedure based on currently accepted techniques, individual judgment and experience. Here is described the placement of the device at a PIP joint. This principle being understood, we can extend its application and principle to other locations. An image intensifier must be made available.

The first Ø1.2 mm pin is introduced transversely in the axis of the epiphysis of P1, ideally at the central point of the axis of rotation of this joint. The second Ø1.2 mm pin is introduced using the drill guide at the diaphysis of P2. By using the FINGERFIX PLIERS forceps, the distal pin will be bent at 90°, on both sides of the finger. This pin will be shortened as shown in the illustrated instructions. Slide the two plastic pins over the proximal pin, one on each side of the finger.

The spring is then "screwed" onto the 90° bent portion of the distal pin as shown in the photos. The spring is advanced halfway. (We can test how the joint reacts to the ligamentotaxis by pulling on the springs). The plastic rod is then introduced into the spring by bending it slightly. We can then complete the "screwing" of the spring to obtain the desired distraction. NB: It is recommended to pull on the spring while tightening it, or while pulling on one of the springs, the other will be tightened and so on. The proximal pin will be twisted 90° and cut on both sides of the finger, making sure to leave enough room for the pin to rotate without touching the skin. Finally, we can cover both ends with protective spindle caps. The spindle protective tips are not mandatory because the FINGERFIX PLIERS give clear and clean cuts. To stabilize the springs, make a frame, using the Ø1.5 mm pin. The frame must be wide enough so that the two ends of the pin can enter one of the two springs while keeping them outside the finger. To maintain the desired distraction, use two pieces of adhesive tape to cover the springs and the U-pin. Finally, indicate on it the number of turns there are between the two pins when you have obtained the optimal distraction.

The amount of distraction to be given will depend on the clinical case and the condition of the joint being treated. We may need a great force on one side and less on the other. Only the surgeon and the use of the image intensifier can best adapt the necessary traction. To avoid dorsal subluxation of the base of P2, a third, more proximal pin can be used.

PRODUCT EXAMINATION AND HANDLING

Prior to use, the product should be visually examined for any signs of contamination or other damage.

CONTAMINATION

ASKORN MEDICAL cannot be held responsible for products contaminated with other substances after the product has left the company.

Surface contamination (talc, dust, dressings, wadding, oils) on the devices can cause a reaction because they are foreign bodies. Handle with care using surgical gloves (rinsed to remove talc) following strictly aseptic technique.

LIMITED WARRANTY, LIMITATION LIABILITY AND WAIVER TO OTHER GUARANTEES

The manufacturer guarantees that all useful precautions have been taken during the manufacture of these products and undertakes to replace any product which, following an investigation by the manufacturer, proves to be defective at the time of delivery. dispatch. Patient selection, surgical procedures, post-surgical treatment and stresses, and device handling are entirely the responsibility of the surgeon. ASKORN MEDICAL has no control over the conditions of use and cannot give any guarantee as to the good or bad effect of the product following the use of the device and cannot be held liable for any loss or damage arising directly or indirectly from the use of this device. The destruction of waste will be done according to the legislation in force.

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

INSTRUCTION FOR USE Of FINGERFIX NON STERILE

For optimal use, the FINGERFIX
PLIERS and the Drill guide are
mandatory




First, put in place proximal 1.2 KW [1]. Then use
disposable drill guide [2] to ensure parallel
placement of 2nd (distal) KW. [3]


Open FINGERFIX PLIERS as wide as possible
and place under the hand.
Bend the distal KW at a 90° angle on both sides of
the finger [4].


Cut bent distal KW and place outwards along the finger. Ee photo above.

Slide plastic rods over proximal KW on both sides of the finger. Select the end of the spring with the longest straight end. Screw each spring completely over both distal KWs and plastic rods as shown [5]. The spring end will lock automatically. If extra traction is required, unscrew the spring away from the plastic rod. Pull distally while re-screwing.

Bend proximal 1.2 mm KW at 90° on each side [6]. Be sure to leave sufficient space to allow plastic rods to rotate without rubbing against the finger. Distribute the space left by FINGERFIX PLIERS on each side equally, as shown on photos [7], [8], [9].

Cut proximal KW as shown, then place KW protective caps as shown and finally crimp them into place.

To stabilise both springs, construct a "frame" using the 1.5 mm KW supplied [0.059in.] creating a U-shape large enough to maintain the springs away from the finger. If necessary, springs can be shortened before constructing the U-shaped frame.

Two strips of tape can be placed around both the springs and the U-shaped frame to secure selected distraction. The number of spring turns between the proximal and distal KWs can be noted on the tapes.

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