

**AKED**.bzh



## **HK2 PLIERS**


REF. 0227 03 00

---

# HK2 PLIERS

THE HK2 PLIER IS SPECIALLY DESIGNED TO:

- » **BEND** laterally Kirschner Wires Ø: 1.2 - 1.5 and 1.8 mm  
Kirschner Wires are not stressed in their bone section during this step.
- » **CUT NEATLY** Kirschner Wires: 1.2 - 1.5 and 1.8 mm
- » **CRIMP** corresponding HK2 connectors.



**BEND** K-Wires Ø:  
1.2 mm - 1.5 mm - 1.8 mm

**CUT NEATLY** K-Wires Ø:  
1.2 mm - 1.5 mm - 1.8 mm

**CRIMP** corresponding HK2  
connectors:  
1.2 mm - 1.5 mm - 1.8 mm

DISTRIBUÉ PAR

**Aked**  
44 C, rue de Bray  
F 35510 Cesson-Sévigné  
Tél. : +33 (0)2 23 35 53 43  
Fax : +33 (0)9 72 26 96 49  
Mail : [contact@aked.bzh](mailto:contact@aked.bzh)

**AKED**.bzh