


AKED.bzh



METAHUS

SINGLE USE EXTERNAL FIXATOR

HK2 METAHUS

SINGLE USE EXTERNAL FIXATOR

ONE HK2 UNIT:

- » 2 KW of same diameter (1.2 mm or 1.5 mm or 1.8 mm)
- » 1 connector to joint 2 KW

One HK2 box contains 10 HK2 units (non sterile) :
20 KW and 10 connectors (by 2, see photo below).

HK218 and HK218F

- Ø 1.8 mm : Distal Radius - Metacarpal bone
- » HK218 : 2 KW
 - » HK218F : 2 KW threaded
- + 1 connector for 2 KW



HK215 and HK215F

- Ø 1.5 mm : Carpal, Metacarpe - Phalange
- » HK215 : 2 KW
 - » HK215F : 2 KW threaded
- + 1 connector for 2 KW



HK212 and HK212F

- Ø 1.2 mm : Carpal - Phalanx
- » HK212 : 2 KW
 - » HK212F : 2 KW threaded
- + 1 connector for 2 KW



HK2 PLIERS [REF. 0227 03 00] WAS SPECIFIQUELY DESIGNED FOR THE USING OF HK2 CONNECTORS.



1

- 1) BENDING the KW
- 2) CRIMPING the connectors
- 3) CUTTING neatly



2



3




HK2

DISTAL RADIUS FRACTURE

Gérard Hoel, MD

University Hospital of St Denis - Reunion Island



METAHUS

BENNETH'S FRACTURE

Philippe Liverneaux, Pr

University Hospital of Strasbourg - France



METAHUS

5TH METACARPAL BASE FRACTURE

Philippe Liverneaux, Pr

University Hospital of Strasbourg - France




METAHUS

5TH METACARPAL NECK FRACTURE

Philippe Liverneaux, Pr

University Hospital of Strasbourg - France





MALLET FINGER

Stéphanie Gouzou, MD
University Hospital of Strasbourg - France

DISTRIBUÉ PAR

Aked
44 C, rue de Bray
F 35510 Cesson-Sévigné
Tél. : +33 (0)2 23 35 53 43
Fax : +33 (0)9 72 26 96 49
Mail : contact@aked.bzh

AKED.bzh