



FINGERFIX

COMPLEX PROXIMAL INTERPHALANGEAL JOINT FRACTURE

FINGERFIX

COMPLEX PROXIMAL INTERPHALANGEAL JOINT FRACTURE

ADVANTAGES

- » Adjustable distraction!
- » Quick and easy to use
- » Ready to use
- » Single use



PACKAGING

- » Single use

CONTENT:

- » 2 plastic wires
- » 2 KW Ø: 1.2mm + 1 KW Ø: 1.5mm (« U » shape frame)
- » 2 springs with open ends
- » 4 KW tip protections
- » 1 single use drill guide to ensure KW parallel placement



One must use our FingerFix plier to ensure ideal setting.



DISTRIBUÉ PAR

Aked

44 C, rue de Bray
F 35510 Cesson-Sévigné
Tél. : +33 (0)2 23 35 53 43
Fax : +33 (0)9 72 26 96 49
Mail : contact@aked.bzh

AKED.bzh