

**AKED**.bzh



**KIROFIX RANGE**

# KIROFIX RANGE

## MALLEABLE DEVICE FOR MAINTAINING HAND POSITION


Made out of a composite bio-compatible material.  
Radio-lucency without lead weight.  
Will not stain surgical gloves.

### ONE WRIST STRAP: TWO POSSIBLE ADJUSTMENTS

Wrist strap / child positioning  
Wrist strap / adult positioning




Supplied with 5 RINGS



KIROFIX Ballast (see accessories)

### KIROFIX Medium - ref. 0226 03 00

Weight: 0.27 kg  
Warranty: 1 year




## ACCESSORIES

### KIROFIX Ballast - ref. 0226 04 00

It is an independent spare part made out of stainless steel. It may be attached to KIROFIX to weigh it.

Weight: 0.780 kg

Warranty: 5 years



### RING - ref. 0226 05 00

RING will be used as an attachment to any KIROFIX.

Sold separately by box of 10 units.

Warranty: 1 year




RING will replace rubber bands made from surgical gloves.



### KIROFIX Large - ref. 0226 01 00

Weight: 0,5kg


Warranty: 1 year



### KIROFIX Small - ref. 0226 02 00

Weight: 0,18kg

Warranty: 1 year



## KIROFIX MEDIUM & CANAL CARPIEN PAR MINI-OPEN

**Pr Philippe Liverneaux**


CHU Strasbourg — France



## KIROFIX MEDIUM & CANAL CARPIEN PAR ENDOSCOPIE

**Dr Didier Fontes**

Clinique du Sport, Paris — France



**Aked**

44 C, rue de Bray

F 35510 Cesson-Sévigné

Tél. : +33 (0)2 23 35 53 43

Fax : +33 (0)9 72 26 96 49

Mail : [contact@aked.bzh](mailto:contact@aked.bzh)